



## RT 300 POCKET CARD

*Please setup the following before fitting the patient onto the bike:*

1. General and patient-specific login
  - Login: "clinician", password: "1776"
  - Input "Patient ID #" (7 digits) and input "PIN #" (4 digits)
  - Revise the patient PIN # when the tablet has an active internet connection
2. Clip the pulse oximeter onto the patient's earlobe. Check that the patient's HR and SpO<sub>2</sub> are visible on the bike's tablet screen before starting the session.
  - Upon session completion, close the pulse oximeter's clip to power down.

*After setting up the patient on the bike:*

3. Ending cycling
  - Pause cycling: press STOP. Pause has a timespan of 15 minutes. Do not press RESTART because this will cause the session to start over again. Press RESUME to continue the session.
  - Stop cycling early with a cooldown: tap COOLDOWN. The therapy session will conclude with a 30 second cooldown and the data will automatically be stored to the tablet.
  - Stop cycling naturally: the therapy session will conclude at 30 minutes with a 30 second cooldown. The data will be stored to the tablet.
  - Stop cycling immediately: press the emergency stop button. Tap "Patient" then "Close Patient" to save the cycling session data onto the tablet.



## RT 300 POCKET CARD

*Please setup the following before fitting the patient onto the bike:*

1. General and patient-specific login
  - Login: "clinician", password: "1776"
  - Input "Patient ID #" (7 digits) and input "PIN #" (4 digits)
  - Revise the patient PIN # when the tablet has an active internet connection
2. Clip the pulse oximeter onto the patient's earlobe. Check that the patient's HR and SpO<sub>2</sub> are visible on the bike's tablet screen before starting the session.
  - Upon session completion, close the pulse oximeter's clip to power down.

*After setting up the patient on the bike:*

3. Ending cycling
  - Pause cycling: press STOP. Pause has a timespan of 15 minutes. Do not press RESTART because this will cause the session to start over again. Press RESUME to continue the session.
  - Stop cycling early with a cooldown: tap COOLDOWN. The therapy session will conclude with a 30 second cooldown and the data will automatically be stored to the tablet.
  - Stop cycling naturally: the therapy session will conclude at 30 minutes with a 30 second cooldown. The data will be stored to the tablet.
  - Stop cycling immediately: press the emergency stop button. Tap "Patient" then "Close Patient" to save the cycling session data onto the tablet.



## TROUBLESHOOTING

1. Ergometer communications error
  - Release emergency stop button
  - Re-insert each connector: power cord, USB plug
2. Patient ID not valid
  - Re-try entering the patient ID
  - Upload patient ID to the tablet while connected to the internet
3. Stop synchronization
  - Press the 'CANCEL SYNCHRONIZATION' box in place of the pulse oximeter icon

Tablet training information video: <https://youtu.be/yjjqfKG68pk>

### Questions

Related to CYCLE? Contact:

Michelle Kho, PT, PhD cell: 289-237-0194

Alex Molloy, cell: 289 808 6947

Related to the tablet? Contact:

RTI Support, 800-609-9166 x316

[icucycle.com](https://www.icucycle.com)

August 21, 2018

Software version 5.02.901



## TROUBLESHOOTING

1. Ergometer communications error
  - Release emergency stop button
  - Re-insert each connector: power cord, USB plug
2. Patient ID not valid
  - Re-try entering the patient ID
  - Upload patient ID to the tablet while connected to the internet
3. Stop synchronization
  - Press the 'CANCEL SYNCHRONIZATION' box in place of the pulse oximeter icon

Tablet training information video: <https://youtu.be/yjjqfKG68pk>

### Questions

Related to CYCLE? Contact:

Michelle Kho, PT, PhD cell: 289-237-0194

Alex Molloy, cell: 289 808 6947

Related to the tablet? Contact:

RTI Support, 800-609-9166 x316

[icucycle.com](https://www.icucycle.com)

August 21, 2018

Software version 5.02.901