The Physical Function in ICU Test - Scored (PFIT-s): Summary

General Notes:

- This test can be performed once a patient is alert enough to cooperate
- If a patient is able to sit in a chair, this position should be used for all components of the test. If they cannot, testing can be completed sitting over the edge of the bed. If they are unable to sit over the edge of the bed, muscle testing can be done in bed, but the patient will be unable to be assessed on the other components of the test
- Conduct the test components in the order they appear on the strength and function case report forms from top to bottom (i.e. strength tests, then sit to stand assistance, then marching on the spot)
- If a test component is conducted more than once throughout multiple attempts to perform the strength and function assessment, use the highest score achieved for PFIT scoring

Testing Procedures:

Test Component	Equipment	Instructions and Additional Considerations			
Shoulder Flexion Strength (grade)	Chair/Bed	If the patient has limited range of motion, test their strength within the available range			
Knee Extension Strength (grade)	Chair/Bed	• If one side (right or left) is stronger than the other, use the stronger side (highest grade) for PFIT scoring			
Sit to Stand Assistance	Chair/Bed	Ideally patient should have their arms crossed at the wrist and held against their chest. If this is not feasible, the patient may keep their arms at their sides or hold on to the arms of the chair.			
		Provide the least amount of assistance required for the patient to safely transfer to the standing position			
Step Cadence	Stopwatch/Clock	The patient may use a walker if required			
(step/min)		Each time a foot hits the floor, it counts as one step (i.e. stepping with the right foot and then the left foot counts as 2 steps)			
		The patient's foot must completely clear the floor with each step for it to count. If the patient's feet do not clear the floor at each step for 6 steps, the test is over			
		If the patient stops marching for longer than 2 seconds the test is over			
		If the patient can march on the spot for 3 minutes, they automatically receive the highest score (3)			

Step Cadence Sample Calculation:

of steps = 75

Time (min: sec) = 1:20 = 60sec + 20sec = 80 sec

Cadence = $\frac{75 \text{ steps}}{80 \text{ sec}} \times \frac{60 \text{ sec}}{1 \text{ min}} = 56 \text{ steps/min}$

Scoring:

PFIT Scoring	0	1	2	3	TOTAL	
Shoulder strength (grade)*	0, 1 or 2	3	4	5		Total PFIT Score
Knee strength (grade)*	0, 1 or 2	3	4	5		
Sit to Stand Assistance	Unable	Assist x2	Assist x1	No Assist		(0 - 12)
Cadence (steps/min)	Unable	>0 - 49	50 - < 80	80+		

References:

[1] Denehy L, De Morton NA, Skinner EH, Edbrooke L, Haines K, Warrillow S, Berney S. A physical function test for use in the intensive care unit: Validity, responsiveness, and predictive utility of the physical function ICU test (scored). Phys Ther. 2013;93(12):1-e3.

As of June 7, 2018